Tubes Radiogènes à Anode Tournante
Drehanoden - Röntgenröhre
Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.

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Manufactured by Varian Medical Systems
Fabrique par Varian Medical Systems
Hergestellt von Varian Medical Systems
Fabricado por Varian Medical Systems

Specifications subject to change without notice.
Spécifications susceptibles d'être modifiées sans préavis.
Technische Daten ohne Gewähr.
Especificaciones sujetas a cambio sin previo aviso.



3 Ø Constant Potential

Abaques de Expositions CEI 60613
Belastungskurven IEC 60613
Diagramas de Exposición IEC 60613

50/60 Hz

0.6 X 0.9 



1.4 X 1.8 

150/180 Hz

0.6 X 0.9 

1.4 X 1.8 

Note:
Rating charts reflect maximum tube performance. Tube operation is ultimately limited by system software.

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

Le Gaine B-225H

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Encaje de B-225H

Maximum Peak Voltage	150 kV
Anode to Ground	75 kV
Cathode to Ground	75 kV
Maximum X-ray Tube Assembly Heat Content	1.5 MJ (2.0 MHU)
Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat)	3,800 W (5,300 HU/sec)
Maximum Heat Exchanger Dissipation	5,000 W (7,050 HU/sec)
Focal Point Position (Central Ray) Within 1mm (X,Y Direction from the center of radiation port.)	
X-Ray Tube Assembly	
Permanent filtration	1.0 mm Al IEC 60522
Loading Factors for Leakage Radiation	150 kV, 20 mA
Federal Standard High Voltage Cable	72
Ambient Air Temperature Limits for Operation	5°C to 40°C
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to 90%
Weight: Housing	35.5 kg (78 lbs)
Heat Exchange	17 kg (38 lbs)
IEC Classification	Class 1
Safety Devices: Thermal Switch - Normally Closed Contact Opening at 85°C ±3°C	
Flow Switch: Normally Open Contacts Contacts close with adequate oil flow.	
Filament Frequency Limits	50 HZ - 20 kHz
Power Supply	DC

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

Note:

1. Heat inputs into housing include tube power, filament power, and stator power.
2. Heating curves based on no restrictions of natural convection around tube housing assembly.
3. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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To install the Varian housing on Picker Scanners certain trunnion mount parts must be replaced. The replacement parts contained in trunnion assembly kit, catalog #CN-184.

The “adjuster” base plate under Picker trunnion ring assembly is retained for use with the Varian trunnion ring base plate.

B-225H Housing and HE-580 Exchanger Installation

Refer to Figure 1 (CN-184 Trunnion Rings)


1.0 Install Base Plate to Heat Exchanger

1.1 Install mounting plate, Item 7, to heat exchanger using spacers and screws, Items 8 and 9, supplied. See Figure #1 (studs on plate face away from heat exchanger). **Locktight is required on all Heat Exchanger hardware.**

2.0 Preassemble Trunnions

2.1 Assemble bottom half of trunnion rings, Items 2 and 10 to base plate, Item 1. Note Item 10 is mounted to front of base plate. See Figure #1. Use four 5/16” x 3/4” cap head screws (Item 5), supplied. Do not tighten.

ILLUSTRATIVE VIEW - FIG. 1

2.2 Temporarily position trunnion mounting bracket and rings on housing. **Use plate trunnion liners.** Tighten four screws (Item 4) holding trunnion rings to mounting bracket allowing small anode - to - cathode motion between rings and housing.

3.0 Installation of Trunnion Mount to Gantry

3.1 Remove trunnion mount from housing and install on gantry “adjuster base plate” with four original 5/8” bolts.

NOTE: Check clearance between base plate and encoder head read cable. If necessary, the cable can be rerouted under the encoder block.

4.0 Adjusting Plate

4.1 In-out adjusting plate, Item 3. Assemble adjusting screw (Item 12), to adjusting tongue, Item 3, supplied and base plate, Item 1. Fasten adjusting tongue to the adjuster base plate with screws, Item 4, supplied.

5.0 Install B-225H Housing

5.1 Position housing into trunnion rings - anode side toward front of gantry. Check very carefully for proper housing and trunnion liner fit into trunnion ring grooves or a housing tilt will occur. Install upper trunnion rings, Items 2 and 10.

6.0 Install HE-580 Heat Exchanger

6.1 Install HE-580 heat exchanger. Use original nuts and lockwashers, **Locktight is required on all Heat Exchanger hardware.**

NOTE: After installation of B-225H Housing and HE-580 Heat Exchanger, verify balance. An additional counterbalance weight may be required.

1	746891	Z AXIS ADJUSTING BLOCK	13
1	746823	ADJUSTING SCREW ASSY	12
2	221091N	5/16 -18 UNC SOC HD CAP SCR x 3” LG	11
1	746580	TRUNNION RING SET (PICKER) - LEFT SIDE	10
5	749642	SPACER	9
5	050268N	1/4 - 20 x 1 LG FHMS ST STL SCR Phillips	8
1	746494	MOUNTING PLATE	7
2	510022N	5/16 -18 UNC SOC HD CAP SCR x 1” LG	6
4	072117N	5/16 -18 UNC SOC HD CAP SCR x 3/4” LG	5
2	221046N	10 - 32 UNF FLAT HD Bolt x 1/2” LG	4
1	746492	PICKER ADJUSTING TONGUE	3
1	745803	TRUNNION RING SET (PICKER) - RIGHT SIDE	2
1	745802	TRUNNION RING MOUNT	1


Stator - Wiring Diagram
Stator - Schéma de Câblage
Stator - Drahtfarbentabelle
Bovina - Diagramas

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WIRING DIAGRAM - FIG. 2

Stator Type:		
Stator Coil Resistance:		
Black to White	14.0 Ohms ±15%	
Green to White	46.0 Ohms ±15%	
Black to Green	60.0 Ohms ±15%	
Starter Voltage:		
50/60 Hz	Start 220 VAC	Run 70 VAC
150/180 Hz	350 VAC	85 VAC
Time to Full Speed:		
50/60 Hz	0 - 3400 RPM	8.5 sec
150/180 Hz	0 - 9500 RPM	18 sec
X-Ray Tube Assembly:		
GS-2096/B-225H	IEC 60601-2-28	

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7.0 Wiring

7.1 Dress stator and heat exchanger cable to “Inboard Drive Chassis”.

7.1.1 Wiring Connections Housing (Also See Figure #2):

TB301		
Black	B	Phase
White	D	Common
Green	C	Phase Shift
Red	E	Metal Center
Blue	Y	Housing Ground

7.1.2 Wiring Connections Heat Exchanger

White	F	AC Power
Black	J	AC Power
Green	Y	Ground
Blue	L	Flow Switch
Red	Connect in series with the white/black wire from thermal switch on housing. Connect second wire from housing thermal switch to TB301 - Z. See Figure #2.	

NOTE: Wire from TB301-E runs to J515-5 (ground) at grid tank (Big Shot system only). The wire will be a current carrier during an exposure. Do not remove J515-5 from ground.

8.0 Heat Exchanger Fuse

8.1 Heat exchanger fuse - replace 2-amp fuse, F106, in remote electronics cabinet with a 5-amp fuse.

9.0 Starter - The GS-2096 can be run at either 10,000 or 3400

RPM. Varian recommends that the GS-2096 be run at 3400 RPM to extend tube bearing life. The following modifications for Big Shot systems are only required if 3400 RPM operations is desired. These changes will include creating an 8-1/2 second rotor boost time for standard speed anode rotation. The changes will give us the required delay for SS04, Rev. C-2 software.

9.1 General X-Ray Interface (Wire Wrap Version, 173507).

9.1.1 Remove the following wire wraps:

- A7 pin 14 to J63
- C11 pin 12 to J57
- B7 pin 2 to A11 pin 13
- A14 pin 12 to A14 pin 5
- A14 pin 5 to A18 pin 5

9.1.2 Add the following wires:

- A18 pin 4 to A14 pin 5
- A18 pin 5 to A14 pin 12

9.2 General Interface (P.C.B. Version, 173533).

9.2.1 Remove the following:

- R52, Q7, R21 and lift U19 pin 13 out of socket

9.2.2 Add the following wire:

- U7-11 to flying lead U19 pin 13

9.3 High Speed Starter

9.3.1 Remove DC brake from circuit.

- 9.3.1.1 Disconnect plugs 520 and 521 from brake cabinet and install directly from power module to high speed starter.

9.3.2 Replace R46 on rotor logic board with 500K trim pot. (**NOTE:** If replacing a Mega HD tube, R21 should be a 500K pot and can be used.)

9.3.3 Adjustments

Starter Voltage	220 VAC
Run Voltage	70 VAC
Start to Run Delay	8 1/2 sec (R46)

9.4 Software

9.4.1 Configure DESC for “Mega HD” and “Hi Speed”

10.0 New Tube Seasoning Procedure - (Must be performed prior to operating tube above 100 kVp.)

10.1 Use this procedure for new tubes or tubes that have been inactive over seven days.

10.2 X-Ray Control Cabinet

10.2.1 Select 80 kV, 150 mA, 6 seconds.

10.2.2 Select large focal spot.

10.2.3 Make 18 exposures - 20 seconds between exposures.

10.2.4 Allow tube to cool **30 Minutes**.

10.2.5 In 10 kV steps, make exposures to 130 kV, 1 exposure at each step.

11.0 Spellman Generator Modifications (See figure 3)

11.1 Connection of Metal Center Section

- 11.1.1 Install GS-2096 in gantry as outlined in the installation procedure of this brochure.
- 11.1.2 Locate J302 in the left side of the gantry.
- 11.1.3 Make connection to J302 (black/white wire).
- 11.1.4 Run about 150 feet or 16 gage wire from the connection made in 11.1.3 to inside the HV chassis of the Spellman Generator. (left side cabinet) (Schematic 440143-001, 440159-001)
- 11.1.5 Locate the Anode HV Module on the right side inside the HV chassis.
- 11.1.6 Look on the backside of module for a screw stud with a blue wire on it which connects to the LV feedback PCB J1-4. (Schematic 340098-001).
- 11.1.7 Add a spade lug to this end of the wire and connect it to the screw stud. This is the **Anode Shield**.
- 11.1.8 This connection will add the metal center section current back into the current sensing network thus preventing a current imbalance error.
- 11.1.9 Make manual exposure at all mA and kVp levels. Verify correct operation, if kVp imbalance occurs continue to Part 11.2.


11.2 kVp Imbalance Adjustment (do only if kVp imbalance occurs).

- 11.2.1 Locate R5, Cathode KVP TP CAL, on the Analog Control BD in the Control Chassis of the Spellman Generator. (Schematic 440147-002 Sheet 2)
- 11.2.2 Adjust R5 until kVp imbalance error is corrected.
- 11.2.3 Make manual exposure at all kVp and mA levels to verify correct operation.
- 11.2.4 If kVp imbalance still occurs, additional turns may need to be made to R5.


11.3 180 HZ Noise Measurement

- 11.3.1 Using the DPS program measure the 180HZ noise.
- 11.3.2 If the 180 HZ noise is out of specification a 68 uf capacitor should be added in parallel with C5 on the Low Voltage Feed Back P.C.B. The negative lead should be connected to E1.

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Note:
1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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