Note: Document originally drafted in the English language.

Product Description

The B-146H housing is designed for Varian rotating Anode inserts having 100 mm diameter targets.

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Manufactured by Varian Medical Systems
Fabrique par Varian Medical Systems
Hergestellt von Varian Medical Systems
Fabricado por Varian Medical Systems


Specifications subject to change without notice.
Spécifications susceptibles d'être modifiées sans préavis.
Technische Daten ohne Gewähr.
Especificaciones sujetas a cambio sin previo aviso.

Product Description	
Nominal X-ray Tube Voltage	125 kV
Maximum Cathode to Ground	63 kV
Maximum Anode to Ground	63 kV
Maximum X-ray Tube Assembly Heat Content	1,750 kJ (2,450 kHU)
Maximum Continuous Heat Dissipation	1000 Watts
2 Liters/minute, 20°C to 40°C fluid	
Maximum Housing Temperature	78°C
X-Ray Tube Assembly	
Permanent Filtration	0.7 mm Al/75 kV IEC 60522
Loading Factors for Leakage Radiation	
125 kV, 8.0 mA	
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Thermal Switch	
Normally Closed	
Open	85°C ±3.9°C (185°F ±7°F)
Contact Rating - 7A @ 120Vac or 30 Vdc max.	

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

Copyright © 2010, Varian Medical Systems. All Rights Reserved.

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Stator Power:
Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. The stator is rated for regular speed and high speed starters. Time to full speed of 102 mm (4 inch) anode series tubes is between 1.3 and 2.0 seconds.

Immediately following high speed anode rotation, the rotor speed must be reduced to 3600 r/min or less within 10 seconds using a suitable dynamic braking device.

No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms.

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