Tubes Radiogènes à Anode Tournante
Drehanoden - Röntgenröhre
Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.

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Note:
1. Limits are based on maximum track rating except for the following codes:
a - Limited by available heat storage.
b - Limited by window heating.
c - Limited by filament emission.
2. H.S. = Heat Storage
kV = Tube Voltage

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Note:
Rating charts reflect maximum tube performance. Tube operation is ultimately limited by system software.

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
Le Gaine B-172H/MX115CT

Das B-172H/MX115CT Gehäuse

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Maximum Peak Voltage	140 kV
Anode to Ground	70 kV
Cathode to Ground	70 kV
Maximum X-ray Tube Assembly Heat Content	2.2 MJ (3.0 MHU)
Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat)	4.0 kW (340 kHU/min)
X-Ray Tube Assembly - IEC 60522	
Permanent filtration	1.2 mm Al at 140kV
Collimator	1.5 mm Al at 140kV
Total Filtration	2.7 mm Al at 140kV
Loading Factors for Leakage Radiation	140 kV, 21 mA
Federal Standard High Voltage Cable	72
Ambient Air Temperature Limits for Operation	5°C to 40°C
Temperature Limits for Storage and Transport	-35°C to +80°C
Humidity	+10% to +90%
Weight: Housing and Heat Exchanger	52.6 kg (116 lbs)
IEC Classification	Class I
Safety Devices - Thermal Switch	
Normally Closed Contact	Opening at 74°C to 79°C
Pressure Switch	
Normally Closed contact	Opens at 5PSIG ± 3PSIG
Power Supply	DC

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

Note:

- Heat inputs into housing include tube power, filament power, and stator power.
- Heating curves based on no restrictions of natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

Stator Type: GE or SARC		
Stator Coil Resistance:		
Black to White	23.5 Ohms ±15%	
Green to White	23.5 Ohms ±15%	
Starter Voltage:	Start	Run
150/180 Hz	385 VAC	85 VAC
Time to Full Speed:		
150/180 Hz	0 - 7500 RPM	4.4 Sec.
X-Ray Tube Assembly:		
GS-1562/B-172H		
GS-1562/MX115CT IEC 60601-2-28		
(Test based on "Maximum Continuous Anode Input Rate")		

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Note:
1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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