Tubes Radiogènes à Anode Tournante
Drehanoden - Röntgenröhre
Tubos de Rayos-X con Ánodo Giratorio

GS-3576HS/MX165CT Housing

Note: Document originally drafted in the English language.

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
Manufactured by Varian Medical Systems
Fabrique par Varian Medical Systems
Hergestellt von Varian Medical Systems
Fabricado por Varian Medical Systems


Specifications subject to change without notice.
Spécifications susceptibles d'être modifiées sans préavis.
Technische Daten ohne Gewähr.
Especificaciones sujetas a cambio sin previo aviso.

Caractéristiques d'Émission du Filament CEI 60613
Kathoden - Emissionskennlinien IEC 60613
Características de Emisión del Catodo IEC 60613

THREE PHASE EMISSION ($\pm .15$ A)
GS-3576HS 0.9 x 0.7

THREE PHASE EMISSION ($\pm .15$ A)
GS-3576HS 1.2 x 1.2

Le Gaine MX165CT

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Encaje de MX165CT

| | |
|--|---------------------------|
| Maximum Peak Voltage | 140 kV |
| Anode to Ground | 70 kV |
| Cathode to Ground | 70 kV |
| Maximum X-ray Tube Assembly Heat Content | 2.6 MJ (3.7 MHU) |
| Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat) | 3.7 kW (5.0 kHU/sec) |
| Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.) | |
| X-Ray Tube Assembly Permanent filtration | 1.2 mm Al IEC 60522 |
| Loading Factors for Leakage Radiation | 140 kV, 26 mA |
| Ambient Air Temperature Limits for Operation | 5°C to 35°C |
| Temperature Limits for Storage and Transport | -34°C to +60°C |
| Humidity | 20% - 70% |
| Weight - Tube and Housing | 75.6 kg (167 lbs) |
| IEC Classification | Class 1 |
| Safety Devices - Thermal Switch, Normally Closed Contact | Opening at 74°C to 79°C |
| Pressure Switch Normally Closed Contact | Opening at 5 PSIG ±3 PSIG |

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
GS-3576HS
TUBE HOUSING ASSEMBLY HEATING AND COOLING

Note:

- Heat inputs into housing include tube power, filament power, and stator power.
- Heating curves based on no restrictions of natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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


| | | |
|----------------------------------|--------------|------------|
| Stator: | | |
| Stator Coil Resistance: | | |
| Black to White | 23 Ohms ±15% | |
| Green to White | 23 Ohms ±15% | |
| Black to Green | 46 Ohms ±15% | |
| Starter Voltage: | | |
| | Start | Run |
| 50/60 Hz | 385 VAC | 75 VAC |
| 150/180 Hz | 450 VAC | 85 VAC |
| Time to Full Speed: | | |
| 50/60 Hz | 0 - 2700 RPM | 10 Sec. |
| X-Ray Tube Assembly: | | |
| GS-3576HS/MX165CT IEC 60601-2-28 | | |

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
ANODE HEATING AND COOLING CURVES

Note:

1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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