

A Review of Intensity Modulated Radiation Therapy (IMRT): A Cost Effective, Personalized Form of Radiation Therapy

Radiation therapy (RT) is used to treat approximately 67% of cancer patients in the US, one or more times during the course of their treatment, or about 1 million patients per year (ASTRO fact sheet, 2009).

In the early 1990s, the standard of care was to deliver RT using a set of intersecting beams, each shaped in two dimensions to the contour of the cancer. Called 3-D conformal radiation therapy (3-D CRT), the combination of two dimensionally shaped beams produces a three dimensional high dose region that approximates the shape of the cancer.

In the late 1990s, the use of a more sculpted beam of RT became common: Intensity Modulated Radiation Therapy, or IMRT. Each IMRT beam does more than simply conform to the shape of the cancer in two dimensions. It actually varies the intensity of the RT according to the dimensions of the cancer in three dimensions. A single IMRT beam is composed of many small “beamlets”; each can be a different intensity. Very complicated cancer shapes can be created, with rapid fall-off of dose immediately outside the cancer. IMRT improves the dosimetry of 3-D CRT, and is often compared with 3-D CRT in clinical studies.

Due to the increase in the utilization of IMRT over the last decade, the complexity of IMRT to plan and administer, and because at first, there was limited published clinical data, some have criticized IMRT as expensive medicine with uncertain benefit. This handout provides an overview of recent data related to the clinical effectiveness of IMRT and describes studies that are currently underway so the reader will understand the reasons for increased adoption of IMRT.

The goal of IMRT is to escalate the dose delivered to the target in order to achieve higher local cancer control rates, without a corresponding increase in normal tissue toxicity. The results are sufficiently encouraging to stimulate further clinical development (see Tables 1 & 2). Initial data indicate that cancer control after IMRT is superior to results from 3-D CRT and randomized controlled phase III trials are underway. The beam shaping and variation in dose intensity of IMRT is done specifically for each cancer patient and is “personalized” RT.

CLINICAL RESULTS SINCE THE INTRODUCTION OF IMRT

This handout reports comparisons between 3-D CRT and IMRT in terms of toxicity, quality of life (QOL) and cancer control outcomes for several different cancers, based on both published data and on clinical trial data not yet published. (It takes an average of 6.75 years for a data set to go from treatment to publication—see Figure 1). For each data set described, the dates of the patient treatments are given. Data have not yet been reported for patients treated after 2006. This handout covers the types of cancer treatments for which the most IMRT data exist. There is no study or data set where patients treated with IMRT have a worse outcome than 3-D CRT with respect to tumor control as measured by local control, disease free survival or overall survival. Not reviewed are data for brain, gastrointestinal, gynecological, sarcoma, and other cancer sites for which there are limited studies.

Extensive data are published on the dosimetric and physics advantages of IMRT compared to 3-D CRT (Braaksmas, 2003; Luxton, 2004).

Breast Cancer

RT is well established as a treatment for early breast cancer patients. Studies report that RT results in both increased local control of the cancer and increased survival rates. IMRT has had a major role in limiting acute and chronic toxicity and improving the quality of life for women who receive RT (Yang, 2009; Poortmans, 2007; Table 3).

Head & Neck Cancer – Nasopharynx

Nasopharyngeal cancer is endemic in China and Southeast Asia. The standard of care is a combination of chemotherapy and RT (Lu, 2009). The toxicities of RT and concurrent chemotherapy are often severe, causing delays or dose reductions during chemotherapy, interruptions of RT, and diminished QOL for patients. When toxicities force alterations in the planned therapy, this can lead to decreased cancer control. The use of IMRT has substantially decreased these toxicities, and decreased interruptions in planned therapy.

The principle results of studies that looked at IMRT as compared with earlier forms of RT for nasopharyngeal cancer patients showed decreased normal tissue toxicity and improved local control of cancers (Table 8). These results are not likely to be unique to nasopharyngeal cancer patients because the cancer biology, pathology and staging are similar in other cancer sites.

Head & Neck Cancer – Oropharynx

A phase III multicenter trial (PARSPORT) in the United Kingdom compared IMRT to 3-D CRT in the treatment of pharyngeal cancer patients (Nutting et al. 2009). The percentage of 3-D CRT patients experiencing grade 2 or worse xerostomia was 64%, compared with 41% from the IMRT group – a statistically significant difference.

A retrospective study (Clavel, 2009) compared the toxicity and the efficacy of 3-D CRT and IMRT administered to patients who were also receiving chemotherapy for locally advanced cancer of the oropharynx. The results after three years of follow-up give significantly improved overall survival, disease free survival, and locoregional control of the tumor with IMRT (Table 5).

French physicians did a matched pair analysis of head and neck patients treated with IMRT vs. 3-D CRT (Graff et al. 2007). They studied 67 pairs of patients. Using validated QOL questionnaires, they reported statistically significant improvements in QOL for patients treated with IMRT, including less dry mouth, sticky saliva, mouth pain, jaw pain, and swallowing and eating difficulties. Xerostomia greater than grade 2 occurred in 67% of the 3-D CRT patients and in only 12% of the IMRT patients. There were no differences in cancer control outcomes.

Lung Cancer

Several physics and dosimetry studies have compared 3-D CRT and IMRT treatment plans for the treatment of locally advanced lung cancer (i.e. stages III and IV). The resulting dose

distributions and dose volume histograms show better sparing of normal tissues with IMRT. The IMRT plans, delivered lower doses to the healthy lung, esophagus, heart, and spinal cord (Christian, 2006; Liu, 2003; Grills, 2003).

Several clinical studies report that higher doses of RT delivered to the cancer result in improved local cancer control (Kong, 2005; Rengan, 2004). Since IMRT makes it possible for physicians to deliver higher doses without causing commensurate levels of toxicity in healthy tissues, future studies may show greater treatment efficacy with IMRT. Current clinical trials are designed to evaluate whether IMRT can deliver higher doses while holding toxicity to acceptable levels. The outcomes that these studies are designed to measure include local cancer control, toxicity, and QOL (Table 3, 4, 5, 8).

Prostate Cancer

The studies reported to date comparing 3-D CRT with IMRT in the treatment of prostate cancer are not controlled trials, but are retrospective comparisons of 2 cohorts of patients treated in different years. Some studies use a “matched pair” form of analysis.

All of these studies concern “early stage” patients. The definition of early stage varies, and involves age, tumor stage, PSA value and Gleason score. Because of these differences, the studies can be somewhat difficult to directly compare. Some also use different criteria for evaluating PSA control as an endpoint.

IMRT in the treatment of prostate cancer is used for two clinical aims: reduction of treatment toxicity and improvement in disease free survival (DFS). In the quest for higher rates of DFS, some centers have used IMRT to escalate the dose to the prostate, delivering doses that would produce unacceptable levels of toxicity using 3-D CRT. Other centers choose to stay with lower doses, and use IMRT only to reduce toxicity. Some level of urinary toxicity from RT to the prostatic urethra, which runs through the center of the prostate, is unavoidable.

Several randomized trials demonstrate that higher doses of 3-D CRT produce a better DFS rate (Goldner, 2009; Zelefsky, 2008). The Fox Chase Cancer Center experience with prostate cancer patients shows a dose response for doses from less than 72 Gy to greater than 76 Gy (Pollack, 2004). IMRT makes it possible to deliver doses that are higher (≥ 80 Gy), and there is evidence that these higher doses produce even longer DFS, especially in low and intermediate risk patients.

Patients in the Memorial Sloan-Kettering Cancer Center (MSKCC) study (Cahlon, 2008) attained a higher PSA control rate and a much lower rectal toxicity rate than patients in the MD Anderson Cancer Center (MDACC) study (Kuban, 2008). This is due to MSKCC patients having been treated with IMRT, while MDACC patients were treated with 3D-CRT (Table 7).

RTOG sponsored a phase III trial (0126) comparing different RT doses in the treatment of prostate cancer. The study, which opened in March 2002 and closed in August 2008, involved 1,532 patients. The participating institutions chose whether to use 3-D CRT or IMRT based on technical capability. Patients were randomly assigned to the 70.2 Gy or 79.2 Gy dose groups after the RT treatment modality was selected. About one-third of the patients received IMRT. The study is in the follow up period, and will report on toxicity levels, QOL, PSA failure rates, disease free survival, and overall survival.

The current RTOG trial (0415) for early stage prostate cancer patients compares two other treatment regimens: delivery of 70 Gy delivered over 28 fractions as compared with 73.8 Gy delivered over 41 fractions. Again, the choice of 3-D CRT or IMRT is made by each institution, prior to the randomly assigned dose regimen. The outcomes measured are the same as in the previous trial described above.

These two RTOG trials will allow for a comparison of 3-D CRT and IMRT.

Ongoing Clinical Trials

A search of the NCI clinical trials database for IMRT trials, conducted in August 2009, yielded 169 results; 55 were not pertinent or had been withdrawn. Table 1 gives a listing by cancer type and study phase. The 21 phase III trials, which are ongoing, will yield more data on the value of IMRT in the ensuing years.

Specific RTOG trials using IMRT are listed in Table 2.

Conclusions

A key question regarding IMRT is whether IMRT is cost effective. The cost is higher than other forms of RT, mostly due to increased physician and other staff time. The additional time spent in planning and delivering IMRT is now a necessary part of maintaining accuracy.

With time, more of the processes involved in planning and delivering IMRT treatments will be computer controlled, so staff time will decrease. Manufacturers of RT equipment and software have made significant progress in producing tools that make IMRT easier and faster to plan and deliver. As clinicians perform more IMRT, their speed and proficiency has improved. New, faster forms of IMRT are now available, which speed up delivery, and thereby increase the throughput and reduce the cost per patient.

The effectiveness is initially evidenced by better treatment plans designed to deliver more dose to the cancer and less to surrounding healthy tissues. This allows for less normal tissue toxicity, which maintains the patients QOL. Toxicity prevention is cost effective when compared to the cost of administering treatments for such toxicities.

Better cancer control occurs when higher doses of RT or better combinations of RT and chemotherapy are possible. This is supported by data in the treatment of head & neck, lung and prostate cancer. No data set reports worse outcomes for IMRT patients (Table 8).

The broad use of IMRT will lead to more clinical studies and opportunities to demonstrate lower cost and increased effectiveness. IMRT use started in 1994 and 15 years later the data on IMRT is starting to mature. Within 5 years the number of studies involving IMRT treatments and their outcomes will double. IMRT is not a treatment offered only in academic medical centers, but is used daily in community-based cancer care settings around the world.

Figure 1 – Date Timetable in Months

It takes an average of 6.75 years for a data set to go from treatment to publication

| Treatment | Follow Up | Data Analysis | Time to Publication | Overall Time  |
|----------------|----------------|---------------|---------------------|--|
| 24 – 36 months | 24 – 36 months | 1 – 3 months | 12 – 24 months | 61 – 99 months (average 81 months or 6.75 years) |

Table 1 – Clinical Trials of IMRT from the NCI Clinical Trials Data System as of August 2009

| Cancer Site | Phase I | Phase II | Phase III | Total |
|----------------------------|-----------|-----------|-----------|------------|
| Breast | 2 | 10 | 4 | 16 |
| Head & Neck | 5 | 28 | 4 | 37 |
| Head & Neck Nasopharyngeal | | 5 | 1 | 6 |
| Lung | 3 | 4 | 2 | 9 |
| Prostate | 5 | 13 | 9 | 27 |
| Anal | | 2 | | 2 |
| Rectal | 2 | | | 2 |
| Pancreatic | | 2 | | 2 |
| Glioblastoma | 1 | | | 1 |
| Brain Mets | 4 | 3 | | 7 |
| Cervical | 1 | 1 | | 2 |
| Sarcoma | | 2 | 1 | 3 |
| Total | 23 | 70 | 21 | 114 |

Table 2 – Radiation Therapy Oncology Group (RTOG) Trials Involving IMRT (N=14) as of December 2009

| RTOG # | Phase | Cancer Site | Therapy | # Patients | Date Started – Completed |
|---------------|-------|-------------------------|---|------------|--------------------------|
| 0615 | II | Head & Neck Nasopharynx | Chemo + 3-D CRT or IMRT | 46 | 12/06 – 12/09 |
| 0522 | III | Head & Neck | Accelerated RT or IMRT + chemo +/- C225 | 942 | 11/05 – 3/09 |
| 0920 | III | Head & Neck | Pre op IMRT +/- C225 | 0 / 700 | 11/09 – |
| 0538 CALGB | III | Small Cell Lung | RT + chemo 3 RT schema | | 3/08 – |
| 0126 | III | Prostate | 70.2 vs 79.2 Gy 3-D CRT or IMRT | 1,532 | 3/02 – 8/08 |
| 0415 | III | Prostate | 70 Gy in 28 fxs vs 73.8 Gy in 41 fxs 3-D CRT or IMRT | 1,067 | 4/06 – 11/09 |
| 0521 | III | Prostate | 72-75.6 Gy 3-D CRT or IMRT | 603 | 12/05 – 8/09 |
| 0621 | II | Prostate | Post op RT 3-D CRT or IMRT | 46 / 76 | 4/08 – |
| 0622 | II | Prostate | Post op 3-D CRT or IMRT + Samarium 153 | 9 / 76 | 4/08 – |
| 0529 | II | Anal | Dose painted IMRT & chemo | 63 | 12/06 – 3/08 |
| 0822 | II | Rectum | IMRT & chemo | 75 | 4/08 – 11/09 |
| 0630 | II | Extremity Sarcomas | Pre op RT 3-D CRT or IMRT +/- chemo | 50 / 102 | 3/08 – |
| 0921 | II | Uterus | IMRT + Cisplatin + Bevacizumab | 0 / 34 | 11/09 – |
| 0418 | II | Cervix, Uterus | IMRT +/- chemo Post op | 48 / 58 | 3/06 – 10/08 |

Table 3 – IMRT Clinical Data
Breast Cancer – 4 Studies

| Reference | Institution / Country | Years Patient Treated | # of Patients | Therapy | Outcomes | Comments |
|----------------|-------------------------|-----------------------|---------------|-------------------------|------------------------|---|
| Barnett, 2009 | Cambridge England | 4/06-6/07 | 404 | 2D Wedged Pair vs. IMRT | | Randomized study Improved dose distribution Final results pending |
| | | | 411 | | | |
| Pignol, 2008 | Canada | 7/03-3/05 | 161 | 2D RT or IMRT | Moist desquamation 48% | Correlates with pain and QOL |
| | | | 170 | | 31% | |
| McDonald, 2008 | Emory | 1/99-12/03 | 124 | 3DCRT or IMRT | Skin Toxicity 52% | No differences in tumor control or survival Retrospective comparison |
| | | | 121 | | 39% | |
| Freedman, 2008 | Fox Chase Cancer Center | 2001-2006 | 405 | 2D RT or IMRT | Skin Toxicity 75% | Retrospective comparison |
| | | | 399 | | 52% | |

Table 4 – IMRT Clinical Data
Nasopharynx Cancer – 5 Studies

| Reference | Institution / Country | Years Patient Treated | # of Patients | Therapy | Outcomes | Comments |
|------------|-------------------------|-----------------------|---------------|----------------|--|--|
| Pow, 2006 | University of Hong Kong | 6/00-7/04 | 21 | 2D RT vs. IMRT | Improved salivary function and QOL scores with IMRT | Stage II |
| | | | 46 | | | |
| Kam, 2007 | Hong Kong | 11/01-12/03 | 28 | 2D RT vs. IMRT | Xerostomia grade ≥ 2 at 1 yr 82% | Better saliva flow rates and QOL with IMRT |
| | | | 28 | | 39% | |
| Lee, 2009 | RTOG | 2/03-11/05 | 68 | IMRT + Chemo | 92% 2 yr local control | |
| Tham, 2009 | Singapore | 2002-2005 | 195 | IMRT +/- Chemo | 3% Xerostomia (grade 3) 90% 3 yr local control 94% 3 yr survival | |
| Wong, 2009 | China | 6/04-12/05 | 175 | IMRT +/- Chemo | 94% local control 87% regional control 87% 3 yr survival | |

Table 5 – IMRT Clinical Data

Oropharynx and Other Head and Neck Cancers – 5 Studies

| Reference | Institution / Country | Years Patient Treated | # of Patients | Therapy | Outcomes | Comments | | |
|---------------|---------------------------------|-----------------------|---------------|-----------------------------|-------------------------------|-----------|----------------------------|----------|
| Nutting, 2009 | United Kingdom Group (Parsport) | 2003-2007 | 47 | 3DCRT vs. IMRT | Xerostomia grade ≥ 2 64% | | Oropharynx | |
| | | | 47 | | 41% p < .05 | | | |
| Clavel, 2009 | Montreal | 1/00-12/07 | 149 | 3DCRT or IMRT RT + Chemo | 3 yr survival | 3DCRT 76% | Stage III or IV oropharynx | |
| | | | 100 | | 3 yr NED survival | 72% | | IMRT 95% |
| | | | | | 3 yr local control | 85% | | 92% |
| Lee, 2006 | Memorial Cancer Center | 9/98-6/04 | 71 | 3DCRT or IMRT | 2 yr. feeding tube | 3DCRT 21% | | |
| | | | 41 | | Xerostomia grade ≥ 2 | 64% | | IMRT 41% |
| van Rij, 2008 | Netherlands | 1/99-12/03 | 88 | 3DCRT or IMRT | IMRT QOL better scores | | | |
| | | | 75 | | | | | |
| Graff, 2006 | France | 1/01-1/05 | 67 | 3DCRT or IMRT | Xerostomia grade ≥ 2 67% | | QOL matched pair analysis | |
| | | | 67 | | 12% | | | |

Table 6 – IMRT Clinical Data

Lung Cancer – 1 Study

| Reference | Institution | Years Patient Treated | # of Patients | Therapy | Outcomes | Comments |
|-----------|-------------|-----------------------|---------------|----------------|---------------------------|---------------------------------------|
| Yom, 2007 | MD Anderson | 8/02-8/05 | 222 | 3D CRT + Chemo | Pneumonitis at 1 year 32% | Stage III, IMRT had larger lung vols. |
| | | | 68 | IMRT + Chemo | 8% | |

Table 7 – Prostate Cancer Therapy for Early Stage Disease

| Treatment | Author | Center | Year Pub. | # of Patients | Median Follow-Up (Months) | Dose (Gy) # Fx | Prostate BED [a/β=1.5] (Gy) | % PSA Relapse | Late Normal Tissue BED [a/β=3] | % Grade 2 + 3 Chronic Urinary Toxicity | % Grade 2 + 3 Chronic Rectal Toxicity | % Chronic Erectile Dysfunction |
|---------------------------|--------------|--|-----------|---------------|---------------------------|-------------------|-----------------------------|---------------|--------------------------------|--|---------------------------------------|--------------------------------|
| 3D Conformal Radiotherapy | Kuban, 2008 | MD Anderson Cancer Center | 2008 | 151 | 104 | $\frac{78}{39}$ | 182 | 22 | 130 | 7+3 | 19+7 | ND |
| High Dose IMRT | Cahlon, 2008 | Memorial Sloan-Kettering Cancer Center | 2008 | 478 | 53 | $\frac{86.4}{48}$ | 190 | 15 | 138 | 13+3 | 2+1 | 34 |
| Prostatectomy | Walsh, 1994 | Johns Hopkins | 1994 | 995 | 4 | NA | NA | 18 | NA | 8 | 1 | 32 |
| Robotic Prostatectomy | Badani, 2007 | Henry Ford Hospital | 2007 | 2766 | 22 | NA | NA | 16 | NA | ND | ND | 21 |

NA = Not Applicable
 ND = No Data Reported

Table 8 – IMRT Clinical Value

| Cancer Type | Decreased Local Toxicity | Increased Tumor Dose | Increased Local Control | Increased Disease Free Survival |
|-------------------|--------------------------|----------------------|-------------------------|---------------------------------|
| Nasopharynx | Yes | NE | Yes | Yes |
| Oropharynx | Yes | NE | Yes | Yes |
| Other Head & Neck | Yes | NE | Yes | NE |
| Breast | Yes | NE | NE | NE |
| Lung | Yes | Yes | NE | NE |
| Prostate | Yes | Yes | Yes | NE |

NE = Not Evaluable

References

1. Badani KK, Kaul S, Menon M. Evolution of robotic radical prostatectomy: assessment after 2766 procedures. *Cancer*. 2007 110(9):1951-8.
2. Barnett GC, Wilkinson J, Moody AM, et al. A randomized controlled trial of forward-planned radiotherapy (IMRT) for early breast cancer: baseline characteristics and dosimetry results. *Radiotherapy and Oncology*. 2009;92:34-41.
3. Braaksma MM, Wijers OB, van Sörnsen de Koste JR et al. Optimisation of conformal radiation therapy by intensity modulation: cancer of the larynx and salivary gland function. *Radiother Oncol*. 2003 Mar;66(3):291-302.
4. Cahlon O, Zelefsky MJ, Shippy A, et al. Ultra-high dose (86.4 Gy) IMRT for localized prostate cancer: toxicity and biochemical outcomes. *Int J Radiat Oncol Biol Phys*. 2008;71:330-7.
5. Christian JA, Bedford JL, Webb S, et al. Comparison of inverse-planned three-dimensional conformal radiotherapy and intensity-modulated radiotherapy for non-small-cell lung cancer. *Int J Radiat Oncol Biol Phys*. 2007;67:735-41.
6. Clavel S, Nguyen D, Després P, et al. Higher dose per fraction and shorter overall treatment time using intensity-modulated radiation therapy versus conventional radiation therapy with concurrent carboplatin and 5-fluorouracil for locally advanced oropharyngeal carcinoma: a comparison of toxicity and efficacy. *J Clin Oncol* 27:15s, 2009 (suppl; abstr 6038).
7. Freedman GM, Li T, Nicolaou N, et al. Breast intensity-modulated radiation therapy reduces time spent with acute dermatitis for women of all breast sizes during radiation. *Int J Radiat Oncol Biol Phys*. 2009;74:689-94.
8. Goldner G, Dimopoulos J, Kirisits C et al. Moderate dose escalation in three-dimensional conformal localized prostate cancer radiotherapy: single-institutional experience in 398 patients comparing 66 Gy versus 70 Gy versus 74 Gy. *Strahlenther Onkol*. 2009 Jul;185(7):438-45.
9. Graff P, Lapeyere M, Desandes E, et al. Impact of intensity-modulated radiotherapy on health-related quality of life for head and neck cancer patients: matched-pair comparison with conventional radiotherapy. *Int J Radiat Oncol Biol Phys*. 2007;67:1309-17.
10. Grills IS, Yan D, Martínez AA, et al. Potential for reduced toxicity and dose escalation in the treatment of inoperable non-small-cell lung cancer: a comparison of intensity-modulated radiation therapy (IMRT), 3D conformal radiation, and elective nodal irradiation. *Int J Radiat Oncol Biol Phys*. 2003;57:875-90.
11. Kam MK, Leung S, Zee B, et al. Prospective randomized study of intensity-modulated radiotherapy on salivary gland function in early-stage nasopharyngeal carcinoma patients. *J Clin Oncol*. 2007;25:4873-9.
12. Khuntia D, Reddy CA, Mahadevan A, et al. Recurrence-free survival rates after external-beam radiotherapy for patients with clinical T1-T3 prostate carcinoma in the prostate-specific antigen era: what should we expect? *Cancer*. 2004 Mar 15;100(6):1283-92.
13. Kong FM, Ten Haken RK, Schipper MJ, et al. High-dose radiation improved local tumor control and overall survival in patients with inoperable/unresectable non-small-cell lung cancer: long-term results of a radiation dose escalation study. *Int J Radiat Oncol Biol Phys*. 2005;63:324-33.
14. Kuban DA, Tucker SL, Dong L, et al. Long-term results of the M. D. Anderson randomized dose-escalation trial for prostate cancer. *Int J Radiat Oncol Biol Phys*. 2008;70:67-74.
15. Lee NY, De Arruda FF, Puri D, et al. A comparison of intensity-modulated radiation therapy and concomitant boost radiotherapy in the setting of concurrent chemotherapy for locally advanced oropharyngeal carcinoma. *Int J Radiat Oncol Biol Phys*. 2006;66:966-74.
16. Lee NY, Harris J, Garden AS, et al. Intensity-modulated radiation therapy with or without chemotherapy for nasopharyngeal carcinoma: radiation therapy oncology group phase II trial 0225. *J Clin Oncol*. 2009;27:3684-90.
17. Liu HH, Wang X, Dong L, et al. Feasibility of sparing lung and other thoracic structures with intensity-modulated radiotherapy for non-small-cell lung cancer. *Int J Radiat Oncol Biol Phys*. 2004;58:1268-79.
18. Luxton G, Hancock SL, Boyer AL. Dosimetry and radiobiologic model comparison of IMRT and 3D conformal radiotherapy in treatment of carcinoma of the prostate. *Int J Radiat Oncol Biol Phys*. 2004 May 1;59(1):267-84.
19. McDonald MW, Godette KD, Butker EK, et al. Long-term outcomes of IMRT for breast cancer: a single-institution cohort analysis. *Int J Radiat Oncol Biol Phys*. 2008;72:1031-40.
20. Nutting C, A'Hern R, Rogers MS, et al. First results of a phase III multicenter randomized controlled trial of intensity modulated (IMRT) versus conventional radiotherapy (RT) in head and neck cancer. *J Clin Oncol*. 27:18s, 2009 (suppl; abstr LBA6006).
21. Pignol J, Olivetto I, Rakovitch E, et al. A multicenter randomized trial of breast intensity-modulated radiation therapy to reduce acute radiation dermatitis. *J Clin Oncol*. 2008;26:2085-92.
22. Pollack et al. Prostate cancer radiotherapy dose response: an update of the Fox Chase experience. *J Urol*. 2004 Mar;171(3):1132-6
23. Poortmans P. Evidence based radiation oncology: breast cancer. *Radiother Oncol*. 2007;84(1):84-101.
24. Pow EHN, Kwong DLW, McMillan AS, et al. Xerostomia and quality of life after intensity modulated radiotherapy vs. conventional radiotherapy for early-stage nasopharyngeal carcinoma: initial report on a randomized controlled clinical trial. *Int J Radiat Oncol Biol Phys*. 2006;66:981-91.
25. Rengan R, Rosenzweig KE, Venkatraman E, et al. Improved local control with higher doses of radiation in large-volume stage III non-small-cell lung cancer. *Int J Radiat Oncol Biol Phys*. 2004;60(3):741-7.
26. Tham IW, Hee SW, Yeo RM et al. Treatment of nasopharyngeal carcinoma using intensity-modulated radiotherapy—the national cancer centre Singapore experience. *Int J Radiat Oncol Biol Phys*. 2009 Dec 1;75(5):1481-6. Epub 2009 Apr 20.
27. van Rij CM, Oughlane-Heemsbergen WD, Ackerstaff AH, et al. Parotid gland sparing IMRT for head and neck cancer improves xerostomia related quality of life. *Radiation Oncology* 2008;3:41.
28. Walsh PC. Radical prostatectomy: a procedure in evolution. *Semin Oncol*. 1994;21(5):662-71.
29. Wong FCS, Ng AWY, Lee VHF, et al. Whole-field simultaneous integrated-boost intensity-modulated radiotherapy for patients with nasopharyngeal carcinoma. *Int J Radiat Oncol Biol Phys*. In press.
30. Yang PS, Chen CM, Lie MC, et al. Radiotherapy can decrease locoregional recurrence and increase survival in mastectomy patients with T1 to T2 breast cancer and one to three positive nodes with negative estrogen receptor and positive lymphovascular invasion status. *Int J Radiat Oncol Biol Phys*. In press.
31. Yom SS, Liao Z, Lin HH, et al. Initial evaluation of treatment-related pneumonitis in advanced-stage non-small-cell lung cancer patients treated with concurrent chemotherapy and intensity-modulated radiotherapy. *Int J Radiat Oncol Biol Phys*. 2007 May 1;68(1):94-102.
32. Zelefsky MJ, Yamada Y, Fuks Z, et al. Long-term results of conformal radiotherapy for prostate cancer: impact of dose escalation on biochemical tumor control and distant metastases-free survival outcomes. *Int J Radiat Oncol Biol Phys*. 2008 Jul 15;71(4):1028-33.



A partner for **life**

USA Headquarters, California

Varian Medical Systems

Palo Alto, CA

Tel: 650.424.5700 | Tel: 800.544.4636

<http://www.varian.com>