


Dimensions are in inches

Front View



Rear View




Image Acquisition Modes (Current)

High-Sense Fluoro:	1024 (h) x 768 (v) 30 fps Gain-4 0.5 pF
Normal Fluoro:	1024 (h) x 768 (v) 30 fps Gain-2 0.5 pF
High-Sense Full-Resolution:	2048 (h) x 1536 (v) 7.5 fps Gain-2 0.5 pF
Dual Gain High Range:	1024 (h) x 1536 (v) 15 fps Gain-1 0.5/4 pF
High-Dose Full-Resolution:	2048 (h) x 1536 (v) 7.5 fps Gain-2 4 pF

Image Acquisition Modes (Digital Board 2.0 with Dynamic-Gain readout required)

Dynamic Gain	1024 (h) x 768 (v)
High Range:	30 fps Gain-1 0.5 pF

Additional Modes:
Consult Varian Medical Systems, Inc.