

## Daily Schedule - Respiratory Gating 1.6

### Monday

Section Description	Start Time	End Time
NOTE: THIS SCHEDULE IS AN EXAMPLE OF WHAT IS COVERED IN THE RESPIRATORY GATING CLASS. CHECK START DATE WHICH MAY NOT BE A MONDAY ....END OF MESSAGE....	08:30 am	10:00 am
Break.	10:00 am	10:15 am
Lecture: Organ motion reduction - comparison of available methods; respiratory gating - amplitude and phase gated treatment; breathing training.	10:15 am	11:30 am
Lunch	11:30 am	1:00 pm
Lecture/Lab: RPM system - system overview, system configuration, control system, file structure, communication with Clinac, CT and simulator.	1:00 pm	2:45 pm
Break	2:45 pm	3:00 pm
RPM system operation and setup overview of Acceptance Procedure, RPM commissioning, basic operation, set-up of system, parameters and data transfer.	3:00 pm	4:30 pm

### Tuesday

Section Description	Start Time	End Time
Lab: Practical Training - simulator. Customer Acceptance procedure. Coaching session.	08:30 am	10:00 am
Break	10:00 am	10:15 am
Lecture: Practical Training - Accelerator. Customer Acceptance procedure - commissioning tests, gated beam delivery, QC/QA tests of RPM. Course ends.	10:15 am	11:30 am