



From: Planning Department

August 09, 2006

Subject: On Board Imager (OBI) Upgrade Memo

The following information should be reviewed with the Customer that is planning to upgrade their existing Clinac or Silhouette Edition Clinac with an On Board Imager (OBI).

**Do not proceed to download an Installation Data Package until the serial number of the Varian Clinac, to be upgraded, has been confirmed with the Customer.**

---

- ALL Clinac 2100, 21EX, 2300, and 23EX with the following serial numbers and **BELOW** will receive/use dedicated hardware only:

21XX - #2799 and below

23XX - # 499 and below

Download the "**Dedicated Hardware - Clinac On Board Imager (OBI) upgrade**" from this web page.

---

- ALL 2100, 21EX, 2300, 23EX and iX with the following serial numbers and **ABOVE** will receive/use integrated hardware only:

21XX - #2800 and above

23XX - # 500 and above

Silhouette Edition Clinac – ALL models

Download the "**Integrated Hardware – Clinac and Silhouette - On Board Imager (OBI) upgrade**" from this web page.

---

For assistance or questions, contact your Varian Planning Representative from this web page, Architectural Support Contact, or the main Varian Planning office.

Varian Medical Systems

Planning Department

911 Hansen Way, M/S C-165

Palo Alto, CA 94304-1028

Phone (800) 278-2747 or (650) 424-5945 • Fax (650) 424-6252

Email: [planning@varian.com](mailto:planning@varian.com)